


# LEVIFLOW<sup>®</sup> Ultrasonic Technology Single-Use High Precision Flowmeters



## LFS-SU Single-Use Flowmeters

LFS-03SU: 0 – 0.8 l/min LFS-10SU: 0 – 20 l/min

.8 l/min Ll 0 l/min Ll

LFS-06SU: 0 – 8 l/min LFS-15SU: 0 – 50 l/min

LFS-20SU: 0 – 80 l/min

## **Ultraclean Non-Invasive Flow Measurement**




Figure 1: Operating principle of ultrasonic single-use sensor




Figure 2: System configuration with single-channel converter for usage with Levitronix<sup>®</sup> Service Software.




Figure 3: System configuration with multi-channel converter.

### INTRODUCTION

The *LEVIFLOW*<sup>®</sup> single-use flowmeters are designed for noninvasive flow measurements of high purity fluids. *Figure 1* illustrates the operating principle. Two piezo-electric transducers, mounted in the sensor housing, generate and receive an ultrasonic wave. The two waves are processed by a signal converter. The difference of the transit time of both waves is proportional to the velocity of the fluid. The wet materials of the single-use sensors fabricated from biocompatible (FDA, USP-VI, BSE/TSE and Animal free) gamma sterilizable polypropylene (PP).

The standard configuration of the *LEVIFLOW*<sup>®</sup> single-use flowmeters (*Figure 2*) consists of a flow sensor and a converter for processing the sensor signals. Various signals (analog output, digital input/output) are provided and can be configured with a PC software. A RS485 bus allows arrays of multiple flowmeters. The sensor value is shown on a 4-digit display. For debugging, data collection and configuration with a PC the *Levitronix*<sup>®</sup> *Service Software* is available at *Levitronix*<sup>®</sup> together with a USB to RS485 adaptor. A stackable 6-channel converter (see *Figure 3*), with almost the same size as the single-channel converter, is available for high volume applications with reduced cabling effort and space need.

### SYSTEM BENEFITS

- High precision (1%) measurement with high turn-down ratio
- No contamination due to non-invasive flow measurement
- No moving parts -> no particle generation
- Improved bubble robustness due to DSP technology
- Flow control together with Levitronix<sup>®</sup> MagLev Pumps
- Easy integration into OEM equipment
- Easy configurable flow sensor parameters (PC software)
- Integrated and configurable totalizer function
- ATEX/IECEx flow sensor versions available.


### **APPLICATIONS**


- High purity and high precision liquid processes
- Sterile non-invasive flow measurement in Pharmaceutical manufacturing
- Biotech processes
- Flow control in combination with Levitronix<sup>®</sup> MagLev pump systems
- Single-use disposable applications


### **SPECIFICATIONS**

Sensor Type Characteristics	LFS-03SU LFS-03SU-SC1	LFS-06SU LFS-06SU-SC1	LFS-10SU LFS-10SU-SC1	LFS-15SU LFS-15SU-SC1	LFS-20SU LFS-20SU-SC1
Flow Range [lpm]	0 - 0.8	0 – 8	0 – 20	0 – 50	0 - 80
Triclamp Fitting Size	3/8" (ID = 6.4 mm)	3/8" (ID = 6.4 mm)	1/2" (ID = 9.4 mm)	1" (ID = 22.2 mm)	1" (ID = 22.2 mm)
Accuracy of Reading Note: Repeatability < Accuracy/2	LFS-03SU > 35 ml/min: ±1% < 1 ml/min: ±10% LFS-03SU-SC1 > 6 ml/min: ±1% < 6 ml/min: ±0.06 ml/min	LFS-06SU: > 1.7 l/min: ±1% < 1.7 l/min: ±17 ml/min LFS-06SU-SC1: > 0.075 l/min: ±1% < 0.075 l/min: ±0.75 ml/min	LFS-10SU: > 4.7 l/min: ±1% < 4.7 l/min: ±17 ml/min LFS-10SU-SC1: > 0.75 l/min: ±1% < 0.75 l/min: ±7.5 ml/min	LFS-15SU: > 10.6 l/min: ±1% < 10.6 l/min: ±106 ml/min LFS-15SU-SC1: > 2 l/min: ±1% < 2 l/min: ±20 ml/min	LFS-20SU: >18.8 //min: ±1 % <18.8 //min: ±188 ml/min LFS-20SU-SC1: > 3.2 //min: ±1% < 3.2 //min: ±32 ml/min
Wetted Surface [cm <sup>2</sup> ] / Vol. [ml] / Weight [g]	29.5 / 4 / 42	32.2 / 4.8 / 42	53.2 / 12.3 / 61	141.2 / 61.7 / 96	173.5 / 125 / 125
Pressure Drop Coefficient C at 20°C ΔP=C x Q <sup>2</sup> , Q=Flow [lpm], ΔP=Press. Drop [kPa]	16.8	0.880	0.0750	0.0101	0.00350
Fluid Temperature / Ambient Temp.	Normal range: 2 – 60 °C ( 35.6 – 140 °F ) / 0 – 40 °C (32 – 104 °F)				
Maximum Fluid Pressure	0 – 0.5 MPa ( 0 – 5 bar, 0 – 72.5 psi)				
Kinematic Viscosity / Sound Speed	0.3 – 40 mm²/s (0.3 – 40 cSt) / 1000 – 2200 m/s				
Wet Materials / Enclosure Classification	Polypropylene (FDA, USP VI, ADI free), Gamma robust for up to 40 kGy / IP-65 (for connected sensor)				
Cable Jacket / Length / Connector	PVC / Various extension cables available. / Circular type (IP-67), lock-release mounting				

Table 1: Specifications of flow sensors (all data based on calibration with water at 20°C)










Dimensions of LFS-03SU

and LFS-06SU







С

D d

LFS-	-10-SU	-15SU	-20SU
Fitting Size Triclamp (ASME 2009)	1/2"	1"	1"
а	138	164.8	231.2
b	80	100	161.2
С	Ø 12.7	Ø 25.4	Ø 25.4
D	Ø 25	Ø 50.5	Ø 50.5
d	Ø 9.4	Ø 22.2	Ø 22.2
е	44.4	79.4	82.5
f	18.6	23	28
g	24	26	34
h	25	30	35
L	110	160	160

Figure 4: Basic dimensions for LFS-SU sensors





Figure 5: ATEX/IECEx flow sensor configuration and components Note 1: Pictures shows LFS-03/06SU with LKM-1.3 mounting kit. Other sizes are similar. Detailed drawings with IFU can be requested at Levitronix®.

Characteristics	Description or Values
Power Supply / Current / Inrush (Start-Up) Current	24 VDC $\pm$ 10% $/$ 150 mA $/$ Peak of 3.8 A within 210 $\mu s$
Ambient Temp / Humidity Range	0 – 40 °C (32 – 104 °F) / 30 - 85% R.H. (no condensation)
Enclosure Classification and Material	IP-20 (indoor use), ABS
Interfaces (See <i>Figure 6</i> for detailed PIN designation and electrical specification)	<ul> <li>RS485 -&gt; MODBUS protocol -&gt; max. array of 99 channels</li> <li>1x Analog Output 4 - 20mA (0 - 20mA configurable)</li> <li>2x Digital Outputs: Flow Alarm, Measurement Error, Volume Counter Pulse, Volume Counter Alarm, Flow as Frequency or Bubble Detection (default: normally open)</li> <li>1x Digital Input: Volume Counter Reset or Zero Adjust</li> <li>4 Digit display (flow rate, error codes), re-zero button</li> <li>Address potentiometers for RS485 address setting</li> </ul>
Configuration Parameters (Available and configurable with RS485/USB converter and configuration software)	Viscosity, Low Cutoff, Dampening constant (filter) Full scale setting, Linearization (15 points), Alarm Outputs (High and Low Alarm) Volume Counter Settings
Weight / Dimensions	130 g / 123 x 75 x 17.5 mm (see <i>Figure 6</i> for details)
Mounting	DIN rail

Table 2: Specification of converter LFC-1C-PC



PIN	Designation	Specification	
1	DC24V+	24 VDC ± 10% Current: 150 mA	
2	DC24V-	Starting: 4.4 A, 2ms	
3	FG	Field Ground	
4	Analog Out +	4 - 20 mA (0 - 20 mA configurable)	
5	Analog Out -	Load Resistance < 600 Ohm	
6	Digital Out1 +	Max. rating: DC30V, 20mA	
7	Digital Out2 +	(open collector) Various configurable options available	
8	СОМ	depending on firmware	
9	Digital In+	Various configurable options available depending on firmware	
10	Digital In-	No-voltage contact or transistor open collector	
11	RS485 +	RS485 with MODBUS	
12	RS485 -	Protocol	

Figure 6: Dimensions and layout of interfaces of single channel converter LFC-1C-PC

### SPECIFICATIONS

Characteristics	6-Channel Converter Type LFC-6CIO-PC		
Power Supply / Current / In-Rush (Start) Current	24 VDC ± 10% / 270 mA / Peak 4.9 within 210 µs		
Ambient Temp / Humidity Range	0 – 50 °C (32 – 122 °F) / 30 - 85% R.H. (no condensation)		
Enclosure Classification and Material	IP-20 (indoor use), ABS		
Interfaces	<ul> <li>RS485 -&gt; MODBUS protocol -&gt; max. array of 99 ch.</li> <li>Stacking of max. 16 converters -&gt; 5 ms DSP process/time per channel</li> <li>4 Digit display (flow rate, error codes), re-zero button</li> <li>Address potentiometers for RS485 address setting</li> <li>6x Analog Outputs: 4 - 20mA (0 - 20mA configurable)</li> <li>6x Digital Outputs: Flow Alarm, Measurement Error, Volume Counter Pulse, Volume Counter Alarm, Flow as Frequency or Bubble Detection (default: normally open)</li> <li>6x Digital Input: Volume Counter Reset or Zero Adjust</li> </ul>		
Configuration Parameters (Available and configurable with RS485/USB converter and service software)	- Viscosity- Linearization (15 points)- Low Cutoff- Alarm Outputs- Dampening constant (filter)(High and Low Alarm)- Full scale setting- Volume Counter and Volume Counter Alarm Settings		
Weight / Dimensions / Mounting	215 g / 140 x 77.3 x 20.5 mm / DIN rail		

Table 3: Specifications for multi-channel converter LFC-6CIO-PC




Figure 7: Dimensions, mounting and stacking concept of converter LFC-6CIO-PC

### ORDER INFORMATION

Pos.	Part Name	Article #	1% Accuracy Flow Range	Fitting	Note
1a	LFS-03SU-Z LFS-03SU-Z-G25 ' LFS-03SU-Z-SC1 LFS-03SU-Z-SC1-G25 ' LFS-03SU-Z-SC1-Ex	100-30375 100-30399 100-30418 100-30419 100-30514	35 – 800 ml/min 6 – 800 ml/min	Triologo 2/0"	
1b	LFS-06SU-Z LFS-06SU-Z-G25 ' LFS-06SU-Z-SC1 LFS-06SU-Z-SC1-G25 ' LFS-06SU-Z-SC1-Ex	100-30377 100-30400 100-30394 100-30406 100-30515	1.7 – 8 l/min 0.075 – 8 l/min	<ul> <li>Triclamp 3/8"</li> </ul>	The LFS-xxSU-Z-Ex sensors are ATEX/IECEx certified and can be placed in Zone 2 explosive atmosphere. Specific mounting kits (see
1c	LFS-10SU-Z LFS-10SU-Z-G25 ' LFS-10SU-Z-SC1 LFS-10SU-Z-SC1-G25 ' LFS-10SU-Z-SC1-Ex	100-30397 100-30405 100-30408 100-30416 100-30516	4.7 – 20 l/min 0.75 – 20 l/min	Triclamp 1/2"	Table 6 Pos 8) are needed for sensor protection. An ATEX/IECEx certified cable sealing system (see Table 6 Pos. 11) is necessary for the sensor cable to leave the Zone 2 area. The sensors have the following marking: 
1d	LFS-15SU-Z LFS-15SU-Z-G25 ' LFS-15SU-Z-SC1 LFS-15SU-Z-SC1-G25 ' LFS-15SU-Z-SC1-Ex	100-30412 100-30111 100-30431 100-30432 100-30517	10.6 – 50 l/min 2 – 50 l/min	Triclamp 1"	(€ 器 區 II 3D Ex mc tc IIIC T80°C Dc
1e	LFS-20SU.1-Z LFS-20SU.1-Z-G25 <sup>7</sup> LFS-20SU-Z-SC1 LFS-20SU-Z-SC1-G25 <sup>7</sup> LFS-20SU-Z-SC1-Ex	100-30483 100-30484 100-30464 100-30465 100-30518	18.8 – 80 l/min 3.2 – 80 l/min	Triclamp 1"	

 Table 4: Standard flow sensor configurations

 Note 1: Gamma irradiated with dosage > 25 kGy.

Pos.	Part Name	Part #	Description	Note / Interfaces
2	LFC-1C-PC	100-30374	Single Channel Converter	Analog Output (4 – 20 mA), 2x Digital Output, 1x Digital Input, RS485 (MODBUS) Note: EMI ferrite (9) for flow sensor cable included in converter package.
3 (A+B)	LFC-6CIO-PC	100-30446	6-Channel Converter with I/O Interfaces (Digital I/O connector 3a and analog output connector 3b included)	RS485 (MODBUS), 6 analog outputs (4 – 20 mA), 6 digital inputs, 6 digital outputs Order Bus Conn. (8a) and Termination Box (8b) as separate article. Note 1: EMI ferrite (9) for bus connection to be ordered as separate article. When stacking multiple converters every sensor cable needs the same EMI ferrite (9).
4 (A-H)	LFC-1C-PC-SK	100-91072	Converter Starter Kit	Flow converter LFC-1C-PC (A) with Ferrite (B), AC/DC desktop supply (C) with international AC mains inserts, sensor cable LFI-C.1-30 (D), converter connection cable LFI-D.1 (E), RS485/USB cable YN-485I-TR (F), USB stick with Levitronix Service Software and product Literature (G).

Table 5: LEVIFLOW<sup>®</sup> converters

Pos.	Part Name	Part #	Features	Special Feature / Description	
5a	LFI-C.1- <b>10/30/60</b>	190-1030 <b>7/8/9</b>	Cable length: 1/3/6 m	Interconnect cable for connection between sensor and converter, PVC jacket	
5b	LFE-C.2- <b>10/30/60</b>	190-1031 <b>0/1/2</b>	Cable length: <b>1/3/6</b> m	Extension cable with wall-mountable connector for cabinet mounting, PVC jacket. Delivered with protective dust cap on wall-mountable connector side.	
6	YN-485I-TR, USB to RS485 Adaptor-TR Isolated	100-30392	Structure/Design	USB connector (A) with termination resistor and cable with connector pair (B and C) for external RS485 wire connection. Magnetically isolated. Cable length is 2m. Included is a USB space saver cable (D).	
7a 7b 7c 7d	Mounting Bracket LMK-1.2 Mounting Bracket LMK-2.2 Mounting Bracket LMK-3.2 Mounting Bracket LMK-4.2	100-91478 100-91479 100-91480 100-91481	Sensor compatibility	For LFS-03SU and LFS-06SU. For LFS-10SU For LFS-15SU For LFS-20SU	
			Material / Sensor Fixation	Anodized Aluminum / Screw type	
8a 8b 8c 8d	Mounting Bracket LMK-1.3 Mounting Bracket LMK-2.3 Mounting Bracket LMK-3.3 Mounting Bracket LMK-4.3	100-91624 100-91625 100-91626 100-91627	Sensor compatibility	For LFS-03SU and LFS-06SU For LFS-10SU For LFS-1SSU For LFS-20SU	
	Note: For Ex applications.		Material / Sensor Fixation	Stainless Steel / Screw type	
9a	Connector Box for LFC-6CIO-PC	100-30447	COMBICON connector	For wiring RS485 and supply of stacks of LFC-6CIO-PC converter.	
9b	Termination Box for LFC-6C	100-30317		For termination of RS485 bus of LFC-6CIO-PC.	
10	LeviFlow Splitting Ferrite	100-30353	EMI filtering of DC supply	For LFC-6CIO-PC supply and bus needed. On flow sensor cables needed in case of stacking of multiple converters.	
11 (A-E)	ATEX Cable Sealing System ACS-A.2 (Roxtec)	100-91628	Sleeve/Nut (A), Gasket (B) Frame (C)	Stainless Steel and EPDMNote: Lubricant (D), measurement platesRoxylon (EPDM rubber)(E) and IFU are included.	

Table 6: Accessories



Figure 8: LEVIFLOW<sup>®</sup> flow sensors and converter







Figure 9: Accessories




Figure 10: Converter starter kit (see Table 2 Position 4) with components

*Levitronix®* is the world-wide leader in magnetically levitated bearingless motor technology. *Levitronix®* was the first company to introduce bearingless motor technology to the Semiconductor, Medical and Life Science markets. The company is ISO 9001 certified. Production and quality control facilities are located in Switzerland. In addition, *Levitronix®* is committed to bring other highly innovative products like the *LEVIFLOW®* flowmeter series to the market.



#### Headquarter and European Contact

Levitronix GmbH Bändliweg 30 CH-8048 Zurich Switzerland

Phone: +41 44 974 4000 E-Mail: <u>salesEurope@levitronix.com</u>

### US Contact

Levitronix Technologies LLC 10 Speen Street, Suite 102 Framingham, Massachusetts 01701 USA

Phone: +1 508 861 3800 E-Mail: <u>salesUS@levitronix.com</u>

#### Japan Contact

Levitronix Japan K.K. Wing Eight 5floor, 4-16-4 Asakusabashi, Taito-ku Tokyo, 111-0053 Japan

Phone: +81 3 5823 4193 E-Mail: <u>salesJapan@levitronix.com</u>

### Taiwan Contact

Levitronix Taiwan 5F, No. 251, Dong Sec. 1, Guangming 6th Rd., Chu Pei City, Hsin-Chu 302, Taiwan, R.O.C.

Phone: +886 3 657 6209 E-Mail: <u>salesAsia@levitronix.com</u>