

PuraLev® Life Science Integrated Flow Controller Series



PuraLev® iF100SU (Single-Use)

Pump Pressure / Flow: 2.0 bar / 17 l/min
 Single-Use Flow Sensor LFS-06SU: 8 l/min
 Single-Use Flow Sensor LFS-10SU: 20 l/min

Low Shear Design - High Cell Viability





Figure 1: Schematic of the main elements of the MagLev centrifugal pump



LEVIFLOW® High Precision Single-Use Flow Sensors




Figure 2: Integrated flow control system with ultrasonic flow sensors

INTRODUCTION

With the *PuraLev® iF100SU* flow control system *Levitronix®* combines its unique magnetic levitation pump technology with its ultrasonic flow measurement technology. The result is a highly integrated precise flow controller with an integrated pressure source. The centrifugal pump, as a pressure source, has no bearings to wear out or seals to break down and fail. The pump impeller is suspended, contact-free, inside a casing and is driven by the magnetic field of the motor (*Figure 1*).

The pump controller, motor and flow converter are integrated into the driver housing. This reduces cabling and setup effort significantly. Fluid flow rate is precisely controlled by electronically regulating the impeller speed without pulsation. The pump head can be easily inserted and removed with an intuitive bayonet socket.

With the lack of mechanical bearings plus the self-contained pump head design, the risk of contamination is drastically reduced. The absence of narrow gaps between the impeller and pump casing, plus the low-shear pump design allows the gentle pumping of sensitive liquids.

SYSTEM BENEFITS

- High precision, dynamics and turndown ratio.
- No dependency on external pressure source.
- Extremely low particle generation due to the absence of mechanically contacting parts.
- Reduced risk of contamination due to the self-contained design with magnetic bearings and ultrasonic technology.
- Very gentle to sensitive fluids due to low-shear design.
- No narrow gaps and fissures where particles or microorganisms could be entrapped.
- Smooth, continuous flow without pressure pulsation.
- No over-pressure situations (compared to roller pumps).
- Biocompatible and gamma sterilizable single-use parts.
- Dry running capability.
- High flow capability with compact design.
- Very low integration costs as no external controller is needed for flow control.
- Proven pump and ultrasonic flow measurement technology.

APPLICATIONS

- Pumping of shear-sensitive liquids and cells.
- Bioprocessing (for example perfusion).
- Recirculation and transfer applications in bioreactors.
- Filtration.

FLOW CONTROL CONCEPT

Figure 5 illustrates the flow control concept. Flow control, pressure generation and flow measurement is done with one unit. This allows realization of sophisticated flow control algorithms and optimizations to various situations.

There is a linear relationship between flow and speed (Figure 4). The speed is precisely controlled with a high resolution over a wide pump speed range. This allows a flow control with high resolution and high turndown ratio compared to non-linear flow control with valve type flow control concepts. Additionally, the highly dynamic speed controller allows fast flow step responses.

As the speed is monitored and the pressure cannot increase uncontrolled at a given speed, there is no need to protect the hydraulic circuit against over-pressure situations as for example for roller pumps in tube clogging situations.

The versatility of *Levitronix*® flow control systems goes far beyond the capabilities of simple flow controllers. In addition to the flow control function, the *Levitronix*® control firmware comes with several condition monitoring features to monitor the integrity of the fluid circuit. *Levitronix*® flow control systems can generate alarms for preventive filter exchange, no-flow conditions or line clogging. Dynamic Condition Trending (DCT) enables failure prediction and scheduling of preventive maintenance

SYSTEM CONFIGURATION – “STAND-ALONE”

Figure 6 and Figure 10 illustrate a “Plug and Play” stand-alone system model with integrated user panel and buttons to set the flow manually. The driver also contains a PLC interface for remote flow control by analog and digital signals. Various accessories are available like a desktop power supply with relevant power cable and signal cables to connect to the PLC.

SYSTEM CONFIGURATION – “EASYCONNECT”

The “EasyConnect” models (see Figure 7 and Figure 12) with according cable accessories are designed to realize various interface configurations with minimal setup effort.

Two Fieldbus connectors (IN and OUT) allow to setup arrays of multiple flow controllers. Therefore, blending configurations as shown in Figure 8 can be realized. The PLC interface allows not only remote control by analog/digital signals but also connections of external sensors hence enabling for example a precise pressure control or monitoring.

SYSTEM CONFIGURATION – “OEM”

The “OEM” models are designed for a compact integration with one integrated hybrid connector containing all available interface signals (see Figure 9 and Figure 14). Basically, all configurations of the “EasyConnect” models are possible allowing the users with integration capabilities to adapt the cable to their needs.




Figure 3: Pressure/flow curves for aqueous liquids (similar to water)




Figure 4: Simplified block schematics of flow control system




Figure 5: Flow control system models

SYSTEM CONFIGURATIONS




Figure 6: Standard "Stand-Alone" system configuration with main accessories
(See section "Order Information" for details to numbered components and other options)




Figure 7: Standard "EasyConnect" system configuration with main accessories
(See section "Order Information" for details to numbered components and other options)

SYSTEM CONFIGURATIONS




Figure 8: Flowcontrol array configuration with "EasyConnect" models for blending applications
(See section "Order Information" for details to numbered components and other options)





Figure 9: Standard "OEM" system configuration
(See section "Order Information" for details to numbered components and other options)

MODEL DESCRIPTION – STAND-ALONE



| Interface | PIN Name | Description | Standard Designation | Hardware Specification |
|---------------------|----------|------------------------------|-------------------------------|---|
| Power Supply | P+ | + 24 VDC | Supply | Voltage: 24 VDC |
| | P- | Power Input Ground / Earth | | Power: 100 W |
| | NC | Not connected. | -- | -- |
| PLC 6 | Ain | Analog Input (Current Input) | Reference (Set) Flow or Speed | Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation) |
| | Ain_GND | Analog In. GND | -- | Reference for Ain |
| | Dout | Digital Output 1 | Status | Open drain, max. 24V, 100mA Reference ground is GND |
| | GND | Analog Ground | -- | Reference for Dout |
| | Din1 | Digital Input 1 | Enable (Reset) | Galvanic separation with optocoupler 2.2 kΩ input resistance, 5-24V for active input |
| | Din_COM | Com. Digi. Input | -- | Reference for digital input. |
| Flow Sensor | 6 Pins | Flow Sensor | -- | Compatible to LEVIFLOW® flow sensors of the LFSC-D and LFS-SU series. |
| Display and Buttons | -- | Display | Flow and Status Display | -- |
| | -- | Up/Down | Setting Flow | -- |
| | -- | On/Off | Enable/Disable | -- |

Figure 10: Interface specifications of standard "Stand-Alone" model





Figure 11: Basic dimensions and description of standard "Stand-Alone" model

MODEL DESCRIPTION – EASYCONNECT



| Connector | PIN Name | Description | Standard Designation | Hardware Specification |
|--------------|----------------|--------------------------------|--|--|
| Power Supply | P+ | + 24 VDC | Supply | Voltage: 24 VDC |
| | P- | Ground / Earth | | Power: 100 W |
| | NC | Not connected. | -- | -- |
| PLC 12 | Dout1 | Digital Output 1 | Status Pump | Open drain, max. 24V, 100mA Reference ground is GND |
| | Dout2 | Digital Output 2 | Status Flow Sensor | |
| | Din1 | Digital Input 1 | Enable (Reset) | Galvanic separation with optocoupler 2.2 kΩ input resistance, 5-24V for active input |
| | Din2 | Digital Input 2 | Zero Adjust | |
| | Din_COM | Com. Digi. Input | -- | Reference for digital input. |
| | Ain1 | Analog Input 1 (Current Input) | Reference Value (Set Flow/Speed) | Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation) |
| | Ain2 | Analog Input 2 (Voltage Input) | Free Configurable | |
| | Ain_GND | Analog In. GND | -- | Reference for Ain1 and Ain2 |
| | Aout1 | Analog Output (Voltage Output) | Actual Flow | 0 – 10V (no galvanic isolation) GND is reference |
| | GND | Analog Ground | -- | Reference for Aout1, Dout1, Dout2 and Pout |
| Pout | Output +24VDC | Supply Output | For supply of external devices (e.g. sensors). (Current 200mA together with Pout or Fieldbus OUT) | |
| NC | Not connected. | -- | -- | |
| GND | Ground | -- | Reference for Pout. | |
| Fieldbus OUT | Pout | Output +24VDC | Supply Output | For supply of external devices (user panels) (Current 200mA together with Pout or Fieldbus OUT) |
| | RS485+ | RS485 + | Field Bus | Modbus protocol |
| | RS485- | RS485 - | | |
| | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. |
| Fieldbus IN | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. |
| | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. |
| | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. |
| | Internal | Internal Bus | Do not connect | Internal bus needed to connect pumps for serial pumping. |
| Flow Sensor | 6 Pins | Flow Sensor | -- | Compatible to LEVIFLOW® flow sensors of the LFS-SU series. |

Figure 12: Interface specifications of standard "EasyConnect" model





Figure 13: Basic dimensions and description of standard "EasyConnect" model

MODEL DESCRIPTION - OEM



| Connector | Designation | Description | Standard Designation | Hardware Specification |
|-------------|-------------|--------------------------------|----------------------------------|--|
| P+ | | + 24 VDC | | Voltage: 24 VDC |
| P- | | Power Input Ground / Earth | Supply | P- to be connected to earth |
| Ain1 | | Analog Input 1 (Current Input) | Reference Value (Set Flow/Speed) | Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation) |
| Ain2 | | Analog Input 2 (Voltage Input) | Free Configurable | Analog voltage input: 0 – 10V (7.9 kOhm, no galvanic isolation) |
| Ain_GND | | Analog Input Ground | -- | Reference for Ain1 and Ain2 |
| Din1 | | Digital Input 1 | Enable (Reset) | Galvanic separation with optocoupler |
| Din2 | | Digital Input 2 | Zero Adjust | 2.2 kΩ input resistance, 5-24V for active input |
| Din_COM | | Common Digital Input | -- | -- |
| Aout1 | | Analog Output (Voltage Output) | Actual Speed | 0 – 10V (no galvanic isolation) GND is reference |
| Dout1 | | Digital Output 1 | Status Pump | Open drain, max. 24V, 100mA |
| Dout2 | | Digital Output 2 | Status Flow | Reference ground is GND |
| GND | | Analog Ground | -- | Reference for Aout1, Dout1 and Dout2 |
| RS485+ | | RS485 + | Field Bus | Modbus RTU protocol |
| RS485- | | RS485 - | Field Bus | Modbus RTU protocol |
| Shield | | Shielding | Shielding | To be connected to earth (see wire No. 2, P-) |
| Flow Sensor | 6 Pins | Flow Sensor | -- | Compatible to LEVIFLOW® flow sensors of the LFSC-D and LFS-SU series. |

Figure 14: Interface specifications of standard "OEM" model

Note 1: Power supply wire cross-section is 1.5 mm² and for signal wires 0.14 mm²

Note 2: For more detailed description of interfaces consult user manual




Figure 15: Basic dimensions and description of standard "OEM" model

FLOW SENSOR SPECIFICATIONS




Figure 16: Basic dimensions for LFS-06SU-Z and LFS-10SU-Z single-use flow sensors compatible with IFD-100 flow control drivers

| Flow Controller Type | PuraLev® iF100SU with LFS-06SU | PuraLev® iF100SU with LFS-06SU-SC1 | PuraLev® iF100SU with LFS-10SU | PuraLev® iF100SU with LFS-10SU-SC1 |
|---|---|---|---|--|
| Characteristics | | | | |
| Flow Range [lpm] | 0 – 8 | 0 – 8 | 0 – 20 | 0 – 20 |
| Accuracy of Reading (at 20°C fluid temperature) Note: Repeatability < Accuracy/2 | > 1.7 l/min: ±1% < 1.7 l/min: ±17 ml/min | > 0.075 l/min: ±1% < 0.075 l/min: ±0.75 ml/min | > 4.7 l/min: ±1% < 4.7 l/min: ±47 ml/min | > 0.75 l/min: ±1% < 0.75 l/min: ±7.5 ml/min |
| Response Time: Step from 10 – 90% of full scale. | < 1s ¹ | < 1s ¹ | < 1s ¹ | < 1s ¹ |
| Fluid Temperature / Ambient Temperature | Normal range: 10 – 60 °C (50 – 140 °F) | | 0 – 40 °C (32 – 104 °F) | |

Table 1: Specifications of flow controller systems PuraLev® iF100SU with LFS-SU single-use flow sensors.

1: Values for to the specific hydraulic circuit optimized flow control parameters. Standard settings are tuned for general flow control stability and may be higher.

ORDER INFORMATION

| System Name | Article # | Pump Head Socket | Flow Control Driver | Note |
|---------------|-----------|------------------|---------------------|--|
| PLD-IF100SU.1 | 100-91482 | PHS-i100.1 | IFD-100.4-02 | OEM - Driver, one multi-purpose connector, pump head socket |
| PLD-IF100SU.2 | 100-91483 | PHS-i100.1 | IFD-100.5-02 | EasyConnect - Driver with interface connectors, pump head socket. |
| PLD-IF100SU.3 | 100-91484 | PHS-i100.1 | IFD-100.6-02 | Stand-Alone - Driver with integrated user panel, pump head socket. |

Table 2: Standard driver system configurations

| Pos. | Component | Article Name | Article # | Characteristics | Value / Feature |
|------|---|---|-----------|--------------------------|--|
| 1a | Single-Use Pump Heads | DCP-200.2 (Triclamp) | 100-90734 | Impeller / Pump Housing | Polypropylene (FDA, USP Class VI, BSE/TSE/Animal free) |
| 1b | | DCP-200.3 (Barb) | 100-90792 | Housing Sealing | Infrared welding |
| 1c | | DCP-200.2-G25 | 100-91078 | In-/Outlet Fittings | Triclamp ½" or Barb ½" for tubing with typical ID 3/8" |
| 1d | | DCP-200.3-G25 (Gamma Irradiated with Dosage ≥ 25 kGy) | 100-91122 | Max. Flow | 21 liters/min / 5.5 gallons/min |
| | | | | Max. Diff.-Pressure | 2 bar / 29 psi |
| | | | | Max. Viscosity | < 20 cP |
| | | | | Wet Pump Volume/Surface | 24 ml / 150 cm ² with Triclamp and 25 ml / 150 cm ² with Barb fittings |
| | | | | Max. Liquid Temp. | 60°C / 140°F |
| | | | | Applicable Sterilization | Gamma radiation up to 40kGy |
| 2a | Integrated Flow Control Driver ("OEM" Model) | IFD-100.4-02 | 100-10138 | Voltage, Power | 24 VDC ±10%, 100 W |
| | | | | Housing | Epoxy coated Aluminum, PP for bottom lid, IP65 ¹ |
| | | | | Interfaces | PLC, RS485 with Modbus protocol (see Figure 14 for details) and flow sensor |
| | | | | Standard Firmware | J3.48 |
| 2b | Integrated Flow Control ("EasyConnect" Model) | IFD-100.5-02 | 100-10139 | Housing | Epoxy coated Aluminum, PP for bottom lid, IP65 |
| | | | | Interfaces | 2x Fieldbus RS485 with Modbus protocol, PLC, power supply and flow sensor |
| | | | | Standard Firmware | J3.48 |
| 2c | Integrated Flow Control ("Stand-Alone" Model) | IFD-100.6-02 | 100-10140 | Housing | Epoxy coated Aluminum, PP for bottom lid, IP65 |
| | | | | Interfaces | User panel with 3 user buttons, PLC, power supply and flow sensor |
| | | | | Standard Firmware | J3.48 |
| 3 | Pump Head Socket | PHS-i100.1 | 100-91053 | Mounting Type | Bayonet type with locking pin |
| | | | | Material | Anodized Aluminum |
| | | | | Assembly Screws | 4 pcs M3 x 8 mm (Stainless Steel, INOX A4) |

Table 3: Specification of standard components
1: Designed and tested for IP67.

| Pos. | Component | Article Name | Article # | Fitting | Wet Material | Note |
|------|-----------------------------------|-----------------------------|-----------|---------------|--|--|
| 4a | LEVIFLOW® Single-Use Flow Sensors | LFS-06SU-Z (8 lpm) | 100-30377 | Triclamp 3/8" | Polypropylene (FDA, USP Class VI, BSE/TSE/Animal free) | See Levitronix® technical brochure of LFS-SU single-use sensor series for more detailed specifications and for other configurations. |
| 4b | | LFS-06SU-Z-SC1 ¹ | 100-30394 | Triclamp 3/8" | | |
| 5a | | LFS-10SU-Z (20 lpm) | 100-30397 | Triclamp 1/2" | | |
| 5b | | LFS-10SU-Z-SC1 ¹ | 100-30408 | Triclamp 1/2" | | |

Table 4: Specification of LEVIFLOW® single-use high-precision (1% accuracy of reading) flow sensors compatible with IFD-100 drivers
Note 1: Extended calibration for wider 1% accuracy range. Note 2: All flow sensors available with gamma irradiation (see LEVIFLOW® product literature for more details).

| Pos. | Component | Article Name | Article # | Characteristics | Value / Feature |
|------|---|--|---|--------------------------------|--|
| 6a | AC/DC Power Supply | TSP 180-124 (Traco) | 100-40018 | Voltage Output / Input | 24 VDC with 180 W / 85 – 132 and 187 – 264 VAC (autoselect) |
| | | | | Basic Dimensions | 110 x 110 x 54 mm (mountable on DIN rail 35 mm) |
| | | | | Certification or Standards | UL/cUL, CB, Semi F47 |
| 6b | Desktop AC/DC Power Supply | GST160A24-R7B IC915 | 100-40020 | Voltage Output / Input | 24VDC, 160W / 85 – 264 VAC, 47-63 Hz |
| | | | | Basic Dimensions | 175 x 72 x 35 mm |
| | | | | Approvals and Country | UL60950-1, CSA C22.2, TUV EN60950-1 |
| | | | | Safety Approvals | Connector for direct connection to power supply of driver with cable length 1.2 m. |
| | | | | Note | |
| 6c | AC Mains Cables (for Desktop power supply 5b) | AMC-1.1 (2m) AMC-1.2 (2.5m) AMC-1.3 (2.5m) AMC-1.4 (2.5m) AMC-1.5 (2.5m) | 190-10331 190-10332 190-10333 190-10334 190-10335 | Approvals and Country | UL, cUL, US, Canada |
| | | | | Approvals and Country | CB, Germany, Denmark, Norway, Finland, Belgium, Netherland, Sweden, Austria |
| | | | | Approvals and Country | PSE, Japan |
| | | | | Approvals and Country | Switzerland |
| | | | | Approvals and Country | CE, United Kingdom |
| 7 | USB to RS485 Adaptor-TR Isolated | YN-485I-TR | 100-30392 | Structure/Design | USB connector (A) with termination resistor and cable (2m) with connector pair (B and C) for external RS485 wire connection. Magnetically isolated. Cable length is 2m. Included is a USB space saver cable (D). Communication over fieldbus of driver with PC. |
| | | | | Purpose | |
| 8 | IPS Cable Power 2 Wires | ICP-2.1-50 (5 m) | 190-10370 | Cable Material / Wires | PVC jacket / 2x 1.5 mm ² |
| | | | | Connection In / Connection Out | Open wires / Circular Intercontec type to driver |
| | | | | Main Purpose | Connection of power supply to "Stand-Alone" and "EasyConnect" drivers |
| 9a | IPS Cable Signal 6 Wires | ICS-1.1-01 (0.1 m) ICS-1.1-10 (1 m) ICS-1.1-30 (3 m) | 190-10343 190-10344 190-10345 | Cable Material / Wires | PVC jacket / 6x 0.08 mm ² and shielding |
| | | | | Connection In / Connection Out | Circular Hirose type / Circular Hirose type |
| | | | | Main Purpose | Fieldbus connection between "EasyConnect" drivers and flow sensor connection. |
| 9b | IPS Cable Signal 6 Wires | ICS-1.2-50 (5 m) | 190-10346 | Cable Material / Wires | PVC jacket / 6x 0.08 mm ² and shielding |
| 9c | | ICS-1.3-50 (5 m) | 190-10389 | Connection In / Connection Out | Connector with screw type plug for open wire connection / Circular Hirose type |
| | | | | Main Purpose | Fieldbus connection to "EasyConnect" driver / To PLC of "Stand-Alone" driver. |
| 10 | IPS Cable Signal 12 Wires | ICS-2.1-50 (5 m) | 190-10347 | Cable Material / Wires | PVC jacket / 12x 0.14 mm ² and shielding |
| | | | | Connection In / Connection Out | Connector with screw type plug for open wire connection / Circular Hirose type |
| | | | | Main Purpose | General connection to PLC of "EasyConnect" drivers. |
| 11 | IPS Cable Hybrid 15 Wires | ICH-1.1-30 (3 m) ICH-1.1-50 (5 m) | 190-10386 190-10341 | Cable Material / Wires | PVC jacket / 2x 1.5 mm ² for supply wired, 13 x 0.14 mm ² for signal and shielding wire |
| | | | | Connection In / Connection Out | Open wires / Circular hybrid connector for driver connection |
| | | | | Main Purpose | General connection integrated driver connector of to "OEM" driver models. |
| 12 | Fieldbus Termination Connector | FTC-1.1 | 190-10348 | Materials | PPS for connector housing and FPM for sealing. |
| | | | | Main Purpose | Termination of fieldbus. |
| 13a | Mounting Kit | LMK-1.2 (for LFS-06SU) | 100-91478 | Material / Structure | Anodized Aluminum / Locking pin concept |
| 13b | | LMK-2.2 (for LFS-10SU) | 100-91479 | Main Purpose | For mounting of LFS-06SU and LFS-10SU flow sensors. |
| 14 | User Panel | LUI-B.1-01 | 100-30448 | Interface / Housing Rating | RS485 / IP65 |
| | | | | Standard Firmware | A3.00 |

Table 5: Specification accessories

ORDER INFORMATION




Figure 17: Flow control system with standard main components



Figure 18: General standard accessories



Figure 19: Standard cables and other accessories

Levitronix® is the world-wide leader in magnetically levitated bearingless motor technology. *Levitronix®* was the first company to introduce bearingless motor technology to the Semiconductor, Medical and Life Science markets. The company is ISO 9001 certified. Production and quality control facilities are located in Switzerland. In addition, *Levitronix®* is committed to bring other highly innovative products like the *LEVIFLOW®* flowmeter series to the market.



Headquarter and European Contact

Levitronix GmbH
Technoparkstr. 1
CH-8005 Zurich
Switzerland

Phone: +41 44 974 4000
E-Mail: salesEurope@levitronix.com

US Contact

Levitronix Technologies Inc.
10 Speen Street, Suite 102
Framingham, Massachusetts 01701
USA

Phone: +1 508 861 3800
E-Mail: salesUS@levitronix.com

Japan Contact

Levitronix Japan K.K.
Wing Eight 5floor, 4-16-4
Asakusabashi, Taito-ku
Tokyo, 111-0053 Japan

Phone: +81 3 5823 4193
E-Mail: salesJapan@levitronix.com

Taiwan Contact

Levitronix Taiwan
5F, No. 251, Dong Sec. 1,
Guangming 6th Rd., Chu Pei City,
Hsin-Chu 302, Taiwan, R.O.C.

Phone: +886 3 657 6209
E-Mail: salesAsia@levitronix.com